

Single Cell RNA-Seq Analysis in Partek® Flow®

HANDS-ON TRAINING


National Institutes of Health
Aug 2019



Xiaowen Wang
Field Application Specialist
Partek Incorporated
support@partek.com

Login and Project Set-up


- Open your preferred web browser (Chrome, Firefox, etc. would work fine)
- Go to the server URL given by your instructor
- Log in using the username and password given to you
- This will open to the Partek Flow homepage
- Click **New Project** and enter project name: SC-RNAseq-[username]
- This will create a new project



Notes: _____

Experiment Description


- 3k peripheral blood mononuclear cells (PBMCs) from a healthy donor
 - Any peripheral blood cell having a round nucleus
- Downloaded from 10X Genomics' dataset repository
 - <https://support.10xgenomics.com/single-cell-gene-expression/datasets/1.1.0/pbmc3k>
- Today, will be importing the filtered gene/cell matrix from this dataset
- *Goal for today: Identify different blood cell populations*
- Partek Flow is versatile, supporting a wide variety of starting file types
- Partek Flow also supports a wide variety of single cell analysis platforms




Notes: _____

Importing Single Cell Data

- Creating a new project automatically opens up the **Data** tab
- To import the data, click **Import data**, then click **Import single cell data**



- Browse to `/home/flow/FlowData/SingleCell_PBMC3k`
- Select all 3 files (2 tsv and 1 mtx), click **Continue**, then click **Next**
- **Note:** Flow also support .h5 output from CellRanger




Notes:

Specify Metadata

- Click the **Use annotation file** checkbox and set the annotation
 - Assembly: *Homo sapiens (human) - hg38*
 - Gene annotation: *Ensembl transcripts release 91*
- Set **Sample name** to *PBMC 2.7K*
- Click **Finish** to import sample. This will create your first data node


<input checked="" type="checkbox"/>	Sample name	Files	Cells	Features
<input checked="" type="checkbox"/>	PBMC3K	Single_Cell_PBMC3k	2700	32643

Annotation

Use annotation file  ☒


Assembly Homo sapiens (human) - hg38 ▼


Gene/feature annotation Ensembl Transcripts release 91 ▼

Feature identifier 


- ☒ Gene (Values: DDX11L1, DDX11L1, WASH7P, MIR6859-1, MIR1302-2HG, ...)
- ☐ Transcript (Values: DDX11L1-202, DDX11L1-201, WASH7P-201, MIR6859-1-20...)
- ☐ gene_id (Values: ENSG00000223972, ENSG00000223972, ENSG00000227232,...)
- ☐ gene_name (Values: DDX11L1, DDX11L1, WASH7P, MIR6859-1, MIR1302-2HG, ...)
- ☐ transcript_id (Values: ENST00000456328, ENST00000450305, ENST00000488147,...)

Counts format

Raw counts  ☒

Report features without counts  ☒

Gene deduplication

Deduplication method 


☒ Mean ☐ Maximum ☐ Sum

Back Finish


Notes:

Analyses Tab Overview and Running a Task

- Go to the **Analyses** tab
- Your first data node, the **Single cell data** node appears
 - *All data nodes are circles*
- Click the data node
- Clicking any node will bring up a **Context sensitive menu** on the right. Only the tasks that can be performed on that node will appear in this menu
- Select **Single Cell QA/QC** from the **QA/QC** section of the task menu



- This runs the **Single Cell QA/QC** task and produces a new task node
 - *All task nodes are rectangles*




Notes: _____


Cell barcode QA/QC

- Check if droplets (barcode) actually contain cells
- Filter out droplets (cell barcode) don't contain cells
- Knee plot: total read count for all the cell, X-axis represent barcodes in decreasing order on total counts

- **Need to be filtered**




- **Good**




Notes: _____

Single Cell QA/AC

- Double click the Single Cell QA/QC task node to open the task report
- *Single Cell QA/QC shows the most popular QC metrics used in the SC genomics community: the number of read counts per cell, detected genes per cell, and % of mitochondrial reads per cell in three violin plots*
- Set the **Read counts filter** to a max of **4000** reads, the **Detected genes filter** to a max of **1200** genes and the **Mitochondrial reads filter** to a max of **5%**
- Click **Apply filter**



- This runs the **Filter cells** task and outputs a new **Single cell data** node



Notes:

Applying a Noise reduction filter

- Click the filtered **Single cell data** node
- Click **Filter features** in the **Filtering** section of the task menu
- This opens the **Filter features** task dialog
- Click the **Noise reduction filter** checkbox
- Create the following filter using the drop-downs and text boxes
 - Exclude features where **value == 0** in at least **99.9%** of the cells
- Click **Finish** to apply the filter

☒ **Noise reduction filter**

Exclude features where value == 0 in at least 99.9 % of the cells

☐ **Statistics based filter**

Filter features by Counts Percentiles

Keep the top 100.0 features with highest variance

☐ **Feature list filter**


Filter features by Include Exclude

i The features in B cells (Values: BLK, CD19, FCRL2, KIAA0125, MS4A1, PNOC...)

Feature identifier i Gene symbol (Values: 5S_rRNA, TSK, A1BG, A1BG-AS1, A1CF, A2M...)

Back Finish


- The **Filter features** task creates a new **Filtered counts** data node



Notes:


Normalizing counts

- Click the **Filtered counts** node
- Click **Normalization** in the **Normalization and scaling** section of the task menu
- Click on the Recommended button
 - **CPM**
 - **Add 1**
 - **Log2**



The screenshot shows a configuration window titled "Normalization order". It features a green "Recommended" button with a thumbs-up icon. Below the title, there is a list of steps: 1. CPM (counts per million), 2. Add 1 (with a text input field containing "1"), and 3. Log 2.0 (with a dropdown menu showing "2.0").

- Click **Finish** to run the **Normalize counts** task



Notes: _____

Identifying Cell Types

- We'll be using a combination of methods to identify some cell types commonly found in PBMCs. Namely:
 - Unbiased clustering (Graph-based)
 - Visualizing expression using
 - Canonical gene markers
 - Gene lists
 - Lassoing cell populations on the t-SNE plot

Cell Type	Gene Markers
T-cells	CD3D, CD3E
Cytotoxic cells	NKG7, GNLY
B cells	CD79A, CD79B (list)
Monocytes	CD68

Notes: _____


Performing graph-based clustering

- Click the **Normalize counts data** node
- Click **Graph-based clustering** in the **Exploratory analysis** section of the task menu
- Click **Finish** to run with default settings



The screenshot shows a configuration panel for the 'Clustering algorithm'. At the top, there is a section titled 'Clustering algorithm' with an information icon. Below it, three radio buttons are visible: 'Louvain' (selected), 'Louvain with refinement', and 'SLM'. Underneath is a section titled 'Advanced options'. In this section, there is a label 'Option set' followed by a dropdown menu currently showing '-- Default --'. To the right of the dropdown is a 'Configure' button. At the bottom of the panel are two buttons: 'Back' and 'Finish'. A mouse cursor is pointing at the 'Finish' button.


- **Graph-based clustering** produces a **Clustering result** data node



Notes: _____

Invoking the t-SNE plot

- Click the **Clustering result** data node
- Click **t-SNE** in the **Exploratory analysis** section of the task menu
- Click **Finish** to run the t-SNE task with default settings
- A **t-SNE** node is produced, double click it to open the t-SNE plot




- We will use the interactive t-SNE plot to view the clustering results and classify our cells
 - *t-distributed stochastic neighbor embedding (t-SNE) is a popular technique for visualizing high-dimensional data*
 - *t-SNE draws cells that are similar to each other across the high-dimensional RNA-Seq data, where each gene is a dimension, close together on the plot*
 - *t-SNE uses principal components analysis to determine which cells are similar to each other*

Notes: _____

Identifying Cell Types


- Select **Graph-based** from the Color by drop-down menu
 - You can see the data has been clustered into 5 clusters
- The **Group biomarkers** table lists genes that distinguish each cluster
- **CD3D**, a T cell marker gene, is listed as a biomarker for Clusters 1 and 5



Notes:

Rotating, panning, and zooming


- To get a better view of the cluster, we can change our view
 - Rotate the plot by left-clicking and dragging
 - Zoom using the mouse wheel
 - Pan by right-clicking and dragging
 - Move the legend by left-clicking and dragging it



Notes: _____

Coloring cells by marker genes


- Click **CD3D** on the *Group biomarkers* table to color by **CD3D** expression
- Cells are now colored by their expression values for **CD3D** from green (0) to red (max)




Notes: _____

Coloring by a second marker gene


- Click the green plus to color by a second gene



- Type **CD3E** in the second text box and select **CD3E** from the list




- The plot is now colored by CD3D (green) and CD3E (red) with cells that express both genes colored yellow



Notes: _____

Viewing expression values for individual cells


- Click a yellow cell on the plot
- The expression values of that cell are listed in the legend
 - Each gene is assigned a color channel (RGB)
 - Cells that express multiple genes have mixed color
- Yellow cells express both CD3D and CD3E



Notes: _____

Selecting by Cluster


- Cells in Clusters 1 and 5 express T cell marker genes, we want to classify these cells as T cells
- Choose **Graph Based** from the **Select by** drop-down menu
- Click the boxes for **1** and **5** to select cells in those clusters
 - *Selected cells on the plot are shown in bold*
 - *The distribution of expression values for selected cells is shown on the legend*




Notes:

Classifying selected cells

- Click **Classify selection**



- Name the classification **T cells**
- Click **Save**
- T cells is added to the **Classifications** section of the menu
 - The number of T cells is listed in parentheses*




- Clear the selection by clicking a blank space on the plot
- Select **Graph-based** from the **Color by** drop-down menu

Notes: _____

Classifying cytotoxic cells


- In the *Selection/Filter* section, choose *Features* tab, type in NKG7 click **+**
- Type expression ≥ 12 to select the cells
- Click Classify selection and name the classification **Cytotoxic cells**
- Save the classification
- Click the plot to clear the selection



Notes:

Coloring by a gene list

- Cluster 3 lists the B-cell marker genes CD79A and CD79B as biomarkers
- To further verify that these are B cells cells, we can use a published list of 92 marker genes for B cells



- Select **List** from the **Color by** drop-down menu


The image is a screenshot of a software interface for visualizing data points. It has three tabs: 'Points', 'Grouping', and 'Layout'. The 'Points' tab is active. Under 'Color by', there is a dropdown menu set to 'List', with a sub-menu showing 'List' and 'B cells'. Below that is a 'Metric' dropdown set to 'PCA'. Under 'Shape by', there is a dropdown set to 'Fixed shape' and a checked 'Filled' checkbox. Under 'Size by', there is a dropdown set to 'Fixed size'. There is a 'Point size' slider set to 4. Under 'Label by', there is a dropdown set to 'Auto'.

- Select **B cells**
- Choose the color metric

Notes: _____

Coloring by a gene list


- Coloring by a list performs principal components analysis on the gene list to identify cells that are distinguished by their expression of genes on the list
- The color of each cell is determined by its value for the first three PCs (PC1 green, PC2 red, PC3 blue)
- The cells from Cluster 3 are colored green and are distinguishable based on their expression of 92 B cell marker genes




- Choose Sum from the Metric drop-down list

Notes: _____

Selecting cells using the 3D lasso tool

- Click the lasso icon to activate the **3D lasso tool** 
- Click and hold to draw a lasso around the cluster of green cells
- Click the starting circle to close the lasso and select the cluster




- Click **Classify selection** and name this group **B cells** (note that some cells may have already been classified)

Notes: _____

Identifying monocytes


- Clear the selection
- Select **Graph Based** from the **Color by** drop-down menu
- Pan and zoom to focus on Cluster 2
- A biomarker for Cluster 2 is a monocyte marker gene, **CD68**, click on this gene in the *Group biomarkers* table



Notes: _____

Classifying monocytes

- Click the lasso icon to activate the **3D lasso tool**
- Click and hold to draw a lasso around the red cells
- Click the starting circle to close the lasso and select the cluster




- Click **Classify selection** and name this group **Monocytes**
- Click the plot to clear the selection

Notes: _____


Viewing classifications


- Reset the view
- Select **Classifications** from the **Color by** drop-down menu
- N/A means the cell doesn't have be classified in the classifications attribute






Notes: _____

Select and filter

- Next we will classify the N/A cells
- In *Selection /Filtering* section *Attributes* tab, choose **Classifications** from the drop-down
- Select **N/A** to highlight the cells in this group
- Click on  to filter only include those cells to view
- Clear filters display all the cells
- Since N/A cells are close to Cytotoxic cells, we will classify them as cytotoxic cells.
- Click **Summary** in the **Classifications** section of the menu to view the **Classifications Summary** table
- Click **Apply classifications** to run the **Classify cells** task



▼ Selection / Filtering 50 selected

   Clear filters

Display options ☒ Dim ☐ Gray

Attributes Features

Select by Classifications ▼

☐ T cells ☐ Cytotoxic cells



☐ B cells ☐ Monocytes

☒ N/A

☐ All

Notes: _____

Exporting visualizations and notebook

- All visualizations in Partek Flow can be saved as publication-quality images
 - To save the t-SNE plot, click the **Save image** button 
- You can also export an image to a **Project Notebook**. The notebook is always associated with a specific project, so specific notes related to the analysis stays with the data
- To send the same image to the notebook, click the **Send to notebook page** button 
- This will prompt you to specify a notebook page to send the image too. Create a new page and call it *Cell classification* and click the **Send** button
- Click the *Cell classification* page link to navigate to the notebook page
- The page has helpful features for writing notes about specific observations, attaching relevant images and can be exported as a pdf file

Notes: _____

Identifying differentially expressed genes

- Now that we have classified cells into cell types, we can compare expression between cell types
- Here, we will compare **Cytotoxic cells** and **T cells** to identify genes that are differentially expressed between these cell types
- Click the **Classified groups** node produced by the **Classify cells** task
- Click **GSA** in the **Differential analysis** section of the task menu
- Click **Classifications** to include it in the GSA model
 - *Adding a factor to the GSA model means that its effects will be considered in the statistical test*
- Click **Next**

Choose which attributes to include in the statistical test

Categorical attributes

Classifications

☒

Graph-based

☐

Numeric attributes

Expressed genes

☐

Mitochondrial reads percent

☐

Total count

☐

Back

Next

Notes: _____

Adding a comparison

- Differential expression analysis lets us compare groups. Here, we want to compare **Cytotoxic cells** to **T cells**
- Click **Cytotoxic cells** in the top panel
- Click **T cells** in the bottom panel
 - *The top panel is the numerator and the bottom panel is the denominator for fold-change calculations*

The image shows a 'Comparison selector' window. It contains two panels, one for the numerator and one for the denominator, separated by 'vs.'. Each panel has a 'Classifications' section with four options: B cells, Cytotoxic cells, Monocytes, and T cells. In the top panel, 'Cytotoxic cells' is selected. In the bottom panel, 'T cells' is selected. At the bottom of the window is an 'Add comparison' button with a mouse cursor pointing to it.


- Click **Add comparison**
- Click **Finish** to run the statistical test
 - Running the GSA task produces a **Feature list** data node

Notes: _____

Viewing results

- Double click the **Feature Lists** data node to open the ANOVA report
- The **Gene list** table in the GSA report lists every gene that was considered by the GSA
 - *Genes are listed starting with the lowest p-value*
- Click the **dots** icon next to **NKG7** under **View** to open a violin plot

	View	Gene ID	Tc
1		NKG7	3.
2		GZMH	
3		GZMB	



- Click **GSA report** to return to the gene list table

Notes: _____

Filtering results

- To identify significantly differentially expressed genes, we can use the **Filter** on the left-hand side of the table
- Set **FDR step up** to **1e-5** and **Fold change** to **-4 to 4**
- The number of genes in the table changes with the filter applied

Filter	
<input type="checkbox"/> Gene ID	◀
<input type="checkbox"/> Total counts	◀
<input type="checkbox"/> P-value	◀
<input checked="" type="checkbox"/> FDR step up	▼
Less than or equ: ▼	1e-5
0	1
<input type="checkbox"/> Ratio	◀
<input checked="" type="checkbox"/> Fold change	▼
From -4	to 4
<input checked="" type="checkbox"/> Exclude range	
<input type="checkbox"/> LSMean	◀
<input type="checkbox"/> Low expressed	◀
<div>Save filter Clear filter</div>	

- Click  to run the **Differential analysis filter** task

Notes: _____

Configuring Hierarchical clustering

- To visualize the differentially expressed genes on our filtered list, we will create a hierarchical clustering heat map
- Click the **Feature list** node generated by **Filter list**
- Click **Hierarchical clustering** in the **Exploratory analysis** section of the task menu
- Uncheck **Cluster samples**
- Check **Filtering** and set to **Include Classification is T cells** OR **Include Classification is Cytotoxic cells**

Cluster samples *i* ☐

Cluster features *i* ☒

☒ **Filtering**

Include samples in selected groups *i*

include ▼ Classifications ▼ is ▼ Cytotoxic cells ▼ OR ✕

include ▼ Classifications ▼ is ▼ T cells ▼ OR ✕

AND

- Under Ordering select **Classifications** from the **Sample order** drop-down menu to order cells by their classification

Ordering

Sample order *i* Classifications ▼

Cytotoxic

T cells ⇅

B cells


Monocytes

- Click **Finish** to run **Hierarchical clustering**


Notes:

Hierarchical clustering heat map

- Double-click the **Hierarchical clustering** node to open the heat map
- Set the **High** value to 2 to balance the colors
- Select **Classifications** from the **Attributes** drop-down menu to label cells with their classification















- Click the save image button to download the heat map as a publication-quality image




Notes:

- We can use **Biological interpretation** tools to learn more about the differentially expressed genes between **Cytotoxic cells** and **T cells**
- Click the **Feature list** node generated after filtering
- Click **Pathway analysis** in the **Biological interpretation** section of the task menu
- Click **Finish** to run enrichment analysis
- Double-click on the **Pathway enrichment** task node to view the report

Gene set 	Description 	Enrichment score 	P-value 	Genes in list 	Genes not in list 	
path.hsa04650	Natural killer cell mediated cytotoxicity	34.85	7.29E-16	19	108	 
path.hsa05332	Graft-versus-host disease	21.15	6.52E-10	9	29	 
path.hsa04060	Cytokine-cytokine receptor interaction	13.64	1.19E-6	15	256	 

- The links on the table open to KEGG pathway maps overlaid with your differential gene expression results



Notes:

Further Training

Self-learning

- Check out <http://www.partek.com/flow-resources> for documentation and additional resources
- Recorded webinars available on <http://www.partek.com/webinars>
- Use the t-SNE to identify additional cell populations in the PBMC 2.7K data. A few suggestions:
 - CD14+ Monocytes
 - CGR3A+ Monocytes
 - CD8A+ Cytotoxic T cells
 - NK cells
- Ready to analyze a multi-sample dataset? Try our Glioma multi-sample tutorial

Regional Technical Support

- Open a support ticket at partek.com/support
- Phone: +1-314-884-6172

Notes: _____
